





Attachment 1: ISC Packing Configurations

RBCs - Short Journey	
	<p>Confirm Upon Receipt:</p> <ul style="list-style-type: none"> → 8L Polyurethane (PUR) Container with lid firmly seated → Two Series 4 Plates. One plate on top of RBCs, the other plate on the bottom → RBCs located below the red line found within the container → Packing Time to the Time of Receipt ≤ 17 hours <p><i>Refer to packing slip for packing time:</i></p> <div data-bbox="836 840 1494 1207"> <p>Packing Slip NUMBER 000000000109  Date and Time Prepared: 03/07/2013 - 10:13 Purchase Order # : SIDE CAMPUS (0000550125) Order Number : 0000000000032 Ordered Qty : 2 Agreed Quantity : 2</p> </div>

RBCs - Long Journey	
	<p>Confirm Upon Receipt:</p> <ul style="list-style-type: none"> → 8L Polyurethane (PUR) Container with lid firmly seated → Six Series 4 Plates to form cube → RBCs located within cube → Packing Time to Time of Receipt ≤ 27 hours <p>Note: If RBC units are packed right after further processing (washed, deglycerolized, phenotype tagging and irradiation), this configuration will be used and the number of RBC units in the ISC will be limited to three.</p>



Platelets



Confirm Upon Receipt:

- ➔ 8L Vacuum Insulated Panel (VIP) Container with vacuum panels intact. See examples below.
- ➔ Six Series 22 Plates to form cube
- ➔ Platelets located within cube
- ➔ Packing Time to Time of Receipt \leq 24 hours

Acceptable Vacuum Panel



➔ Panels are rigid and skin is tight.


Unacceptable Vacuum Panel



➔ The loose skin and lack of rigidity of the back panel indicates the panel has lost its vacuum.




Frozen Products



Confirm Upon Receipt:

- 8L Polyurethane (PUR) Container with lid firmly seated
 - Six Series 20M Plates to form a cube.
 - Packing Time to Time of Receipt ≤ 27 hours
- Plates will be at -20°C or colder at the time of receipt.
Gloves should be used to handle plates.




Confirm Upon Receipt:

- Expanded Polystyrene (EPS) Container
 - Six Series 20M Plates to form a cube. The top plate will be placed directly on the frozen products. Packing material will be placed above the top plate, when required.
 - Packing Time to Time of Receipt ≤ 27 hours
- Plates will be at -20°C or colder at the time of receipt.
Gloves should be used to handle plates.

Note: In contingency situations, dry ice may be used instead of Series 20M Plates.



Room Temperature Plasma Protein Products (PPP) - 16 Litre




Confirm Upon Receipt:

- 16L Polyurethane (PUR) Container with lid firmly seated
- Two square Series 4 Plates. One plate will be above product and one plate below product
- Packing Time to Time of Receipt ≤ 27 hours

Note: Room temperature PPP may also be received in the original boxes if they are transferred in Canadian Blood Services' vehicles where acceptable shipment temperatures can be confirmed.

Room Temperature Plasma Protein Products (PPP) - 8 Litre




Confirm Upon Receipt:

- 8L Polyurethane (PUR) Container with lid firmly seated
- Six Series 4 Plates to form cube
- PPP located within cube
- Packing Time to Time of Receipt ≤ 27 hours




Refrigerated Plasma Protein Products (PPP) - 16 Litre



Confirm Upon Receipt:

- 16L Polyurethane (PUR) Container with lid firmly seated
- Three square Series 4 Plates. One plate will be on the bottom, one plate will be on top of the insert and the third plate at the top of the ISC
- Note:** Product can be placed in the insert between the bottom plate and the middle plate and between the middle and top plate. If product is only packed in the insert, the middle and top plates will be together on top of the insert.
- Packing Time to Time of Receipt ≤ 27 hours.
- RBC units can also be placed in this packing configuration, either separately or with refrigerated PPP.

Refrigerated Plasma Protein Products (PPP) - 8 Litre



Confirm Upon Receipt:

- 8L Polyurethane (PUR) Container with lid firmly seated
- Six Series 4 Plates to form cube
- PPP located within cube
- Packing Time to Time of Receipt ≤ 27 hours
- RBC units can also be placed in this packing configuration, either separately or with refrigerated PPP.