

Brampton National Platelet Immunology Reference Laboratory (NPIRL) Platelet / HLA Testing Services

Neonatal Platelet Investigation – Paternal NPIRL_PL-03

TEST DESCRIPTION

The Fetal / Neonatal Alloimmune Thrombocytopenia (FNAIT) paternal Investigation includes HPA genotyping and platelet crossmatching when paternal platelets are available.

SPECIMEN AND REQUISITION REQUIREMENTS

Specimen(s)

- Five (5) full 5-7 ml EDTA (lavender) tubes mixed thoroughly by gentle agitation.
- Label specimen(s) with the required minimum information: patient's last name, first name, PHN/ (PHIN) or hospital number or other unique identifier, date of collection and facility name.

Complete Requisition (must include)

- · Patient's last name, first name, date of birth and PHN/(PHIN) or hospital number or other unique identifier
- Facility Name
- · Physician/Health Care Provider name
- · Phlebotomist name, classification, initial
- · Date/time of collection
- · Name, facility, address, fax and contact number of individual to whom the report will be sent

Requisition(s)

Platelet Immunology Requisition_NPIRL

PRE-SHIPPING STORAGE

- If unable to ship within 24 hours of collection store refrigerated at 1-10°C. Refrigerated samples must arrive within 48 hours of collection.
- If unable to ship within 48 hours of collection; separate and freeze serum sample only.
- Paternal samples for platelet crossmatching against maternal serum must arrive within 72 hours of collection (preferably 48 hours).
- · Do not open EDTA Samples

SHIPPING INSTRUCTIONS

 Contact the National Platelet Immunology Reference Laboratory (NPIRL) prior to shipping samples so that information may be obtained to assess the priority of the requested investigation.

Shipping

- Ship in a container with an ice pack (gel pack stored at 2 8° C may be used for refrigerated samples).
- · Select shipping method for container to arrive at testing site within 24 hours.
- Fax shipping information (way bill #) to 905-494-8149. Indicate the patient name and type of investigation on the fax.



Brampton National Platelet Immunology Reference Laboratory (NPIRL) Platelet / HLA Testing Services

Neonatal Platelet Investigation – Paternal
NPIRL_PL-03

SEND TO

Canadian Blood Services
Brampton Centre
National Platelet Immunology Reference Laboratory (NPIRL)
100 Parkshore Drive
Brampton, ON L6T 5M1
Tel: 905-494-5257

Fax: 905-494-8149